


# FlashTite® Immediate Use Sterilization Containers

For safe and effective, closed immediate use steam sterilization, use SteriTite® sealed containers with filterless FlashTite® valve plate. Order the FlashTite container with its respective basket. This sealed system meets TJC, AAMI, AORN and FDA guidelines and may be used for gravity displacement and pre-vacuum steam immediate use sterilization.



Part #	Description	Dimensions (L x W x H)
<b>SC06HGV</b>	6" High Half Size FlashTite Sealed Container w/lid	12.3"x11"x6" (31.2x28x15.24 cm)
<b>SC08HGV</b>	8" High Half Size FlashTite Sealed Container w/lid	12.3"x11"x8" (31.2x28x20.3 cm)
<b>SC06QGV</b>	6" High Mid Size FlashTite Sealed Container w/lid	18.1"x11"x6" (46x28x15.24 cm)
<b>SC08QGV</b>	8" High Mid Size FlashTite Sealed Container w/lid	18.1"x11"x8" (46x28x20.3 cm)
<b>SC06FGV</b>	6" High Full Size FlashTite Sealed Container w/lid	23.5"x11"x6" (59.7x28x15.24 cm)
<b>SC08FGV</b>	8" High Full Size FlashTite Sealed Container w/lid	23.5"x11"x8" (59.7x28x20.32 cm)

SteriTite is a rigid, reusable sterilization container system, constructed of a heat treatable grade of anodized aluminum and passivated stainless steel hardware for a long useful life. Ideal for standardization, the SteriTite container is the only sterilization system that is corrosion resistant for all sterilization modalities, and the container of choice for custom designed instrument sets.